

# **Beam Monitoring Data**

**Brett Viren** 

bv@bnl.gov

Physics Department
Brookhaven National Lab

# **Beam Monitor Data Channelling Process**



- Primary proton and secondary hadron/muon monitors
- Reads from "ACNET" via an XML-RPC intermediate server
- Writes RawBlocks to rotorooter like DAQ or DCS

### Glossary



- ACNET: hardware/software reading out FNAL beam devices. VMS based and difficult to use. (Used at BNL ATF beam tests).
- XML-RPC: Protocol for doing Remote
   Procedure Calls by using an XML (text) data
   format on top of HTTP "POST" requests.
   Implementations in Perl (trivial and powerful)
   and C++ (less so)

### Beam Division's XML-RPC System



- Provides 2 URLs:
  - For simple query/response:

http://www-bd.fnal.gov/xmlrpc/Accelerator Post list of devices, get their values in the response.

For registering callbacks:

http://www-bd.fnal.gov/xmlrpc/Export
Post a list of devices, a trigger devices and the URL
of another XML-RPC server to call back with the
results when triggered.

Acts as a proxy for reading from ACNET

#### **Status**



- BeamData package
- Requirements doc (→)
- Test implementation working
- Currently implemented in Perl + SWIG wrappers around libRotoTalk.

### **Operation**



- 1. Start rotorooter
- 2. Start bmchanneld
- 3. bmchanneld connects to Beam Divisions XML-RPC server (BDX) with a list of readout devices, a trigger and a URL.
- 4. bmchanneld goes into listen mode (becomes an XML-RPC server) on that URL
- 5. BDX calls back each time a trigger occurs.
- 6. bmchanneld receives BDX POSTs
- 7. bmchanneld queries for out of spill data.

# Configuration



- List of in-spill and trigger devices
- List of out-of-spill devices
- Delay after trigger to read OOS devices
- Trigger period to read OOS devices
- Hosts, URLs, ports, etc.

#### **Issues**



- Array type devices partially implemented on XML-RPC end. Not tested.
- No offline side classes written
- Only tested with simple devices, actual devices still being definition.
- Not much response to requirements doc.